

# Wranglr instructions

## Sample information file tab

Input sample information in a .xlsx file (see specification file). Wranglr will process the sample information file and produce files that are compatible with MPP and Worklist.

Wranglr provides useful information to check that your sample information file meets the requirements and that the file is read correctly.

Wranglr	Sample information file	MPP files
---------	-------------------------	-----------

Validate sample information file:

Input sample info file **1)**

Browse...

test.xlsx

Upload complete

Warnings: **2)**

- First column is not named "SAMPLE\_ID"
- Make sure you have no special characters like (), [] or spaces in column names. This warning might also be due to duplicated column names.
- Check if a column name starts with a number (or X[number])!

Number of groups: 2   **3)**

Number of time points: 1

Dimensions: 168 (rows) x 7 (cols). Total observations: 1176   **4)**

Column classes:

- factor: 2, numeric: 5

Warnings (amount of columns with possible errors):

- duplicated\_columns: 1

- 1) Input Excel sheet as described in sample\_form\_specification.docx
- 2) Wranglr displays warnings if your file does not meet the requirements. Any warnings here will prevent processing and must be dealt with.
- 3) If column(s) named GROUP and/or TIME are found, the number of groups and/or time points found will be displayed.
- 4) Some extra information on the sample information sheet. Warnings here do not prevent processing

## Randomize and add QC

### Project title:

### Project code

Project code can be chosen freely (Metavuo will be used to generate project codes in the future.) Wranglr will suggest a letter or combination of letters that it hasn't seen before

☐ Save project code

Choose this to save the project code, so Wranglr doesn't suggest it to anyone else. Wranglr also remembers the code in the case you rerun the same project.

### Project home folder

Project home folder, will be concatenated to the beginning of data file paths

### Sample run order

Choose method for determining run order: original order, globally random or random within groups (e.g. sample types)

### Column containing groups for randomization

If "Random within group" is chosen, a selection for the group column will appear

☒ Run samples from same subject sequentially

If multiple samples originate from the same subject, they can be run sequentially to reduce variance between them

### Subject ID column

If the above checkbox is ticked, subject ID column can be chosen

### Sample positioning type

☒ 96-well plate

☐ 54-vial plate

### Number of plates available

Choose how many plates are available to the machine. (the new Orbitrap can take 3 plates, while the older machine can only take 2)

### Number of samples before every QC:

Choose the interval of the QC samples

If you have run the project with the same title before, Wranglr remembers the project code you used, and suggests it as a default.

Running samples from the same subject sequentially does not interfere with the randomization order. For example, if Sample run order is set to "Global Random" and the "Run samples from same subject sequentially"-checkbox is ticked, the subject order will be randomized, but all samples from the same subject will be run sequentially (in random order).

If Sample run order is set to "Random within group" as above and the "Run samples from same subject sequentially"-checkbox is ticked, the subject order is randomized inside the groups and all samples from the same subject will be run sequentially (in random order).

### Modes

☒ HILIC neg ☒ HILIC pos ☒ RP neg ☒ RP pos Choose which modes will be run

### Number of QC before beginning of mode:

HILIC neg

10

HILIC pos

10

RP neg

10

RP pos

10

**QC sample position of mode:** The same position will be given to every QC sample in the respective mode. Defaults to the last position on the second plate.

HILIC neg

P2-H12

HILIC pos

P2-H12

RP neg

P2-H12

RP pos

P2-H12

☒ Add AutoMSMS to Worklist files

☒ Run AutoMSMS on a QC sample

### Choose samples for AutoMSMS

4 7

Choose non-QC samples for AutoMSMS by original sample ID

### Second column for MPP

GROUP

Choose the column used for grouping in MPP. Defaults to GROUP column if found, else defaults to QC.

Submit & Process

Z

Note that the number of QC samples in the beginning of modes can change in course of analysis. In this case, instead of manually changing the worklist file, you can rerun the processing with corrected values. The randomization of run order calculates its seed number from the project title, so the run order will be identical if the sample information file of any individual project is reprocessed.

Example of a processed sample information file used in MPP and further analysis:

Preview:

☒ Processed sample information file for MPP

☐ Worklist file

Show  entries

DATAFILE	GROUP	RUN_ORDER	INTERNAL_SAMPLE_ID	SAMPLE_ID	QC
Demo_HILIC_neg_0001		1	QC_1		TRUE
Demo_HILIC_neg_0002		2	QC_2		TRUE
Demo_HILIC_neg_0003		3	QC_3		TRUE
Demo_HILIC_neg_0004		4	QC_4		TRUE
Demo_HILIC_neg_0005		5	QC_5		TRUE
Demo_HILIC_neg_0006		6	QC_6		TRUE
Demo_HILIC_neg_0007		7	QC_7		TRUE
Demo_HILIC_neg_0008		8	QC_8		TRUE
Demo_HILIC_neg_0009		9	QC_9		TRUE
Demo_HILIC_neg_0010		10	QC_10		TRUE
Demo_HILIC_neg_0011	C	11	A_1	115	FALSE
Demo_HILIC_neg_0012	B	12	A_2	42	FALSE
Demo_HILIC_neg_0013	A	13	A_3	68	FALSE
Demo_HILIC_neg_0014	B	14	A_4	56	FALSE
Demo_HILIC_neg_0015	C	15	A_5	75	FALSE
Demo_HILIC_neg_0016	B	16	A_6	74	FALSE
Demo_HILIC_neg_0017	B	17	A_7	5	FALSE
Demo_HILIC_neg_0018	B	18	A_8	58	FALSE

Wranglr generates four new columns: DATAFILE, RUN\_ORDER, INTERNAL\_SAMPLE\_ID and QC. The internal sample IDs begin with project-specific code. The processed sample information file will also include all the columns of the original sample information file.

## Example of a worklist file:

### Preview:

- ☐ Processed sample information file for MPP  
☒ Worklist file

Show  entries

Search:

Sample Name	Sample Position	Method	Data File
QC_1	P2-H12	HILIC_neg.m	D:\MassHunter\Data\Demo\Demo_HILIC_neg\Demo_HILIC_neg_0001.d
QC_2	P2-H12	HILIC_neg.m	D:\MassHunter\Data\Demo\Demo_HILIC_neg\Demo_HILIC_neg_0002.d
QC_3	P2-H12	HILIC_neg.m	D:\MassHunter\Data\Demo\Demo_HILIC_neg\Demo_HILIC_neg_0003.d
QC_4	P2-H12	HILIC_neg.m	D:\MassHunter\Data\Demo\Demo_HILIC_neg\Demo_HILIC_neg_0004.d
QC_5	P2-H12	HILIC_neg.m	D:\MassHunter\Data\Demo\Demo_HILIC_neg\Demo_HILIC_neg_0005.d
QC_6	P2-H12	HILIC_neg.m	D:\MassHunter\Data\Demo\Demo_HILIC_neg\Demo_HILIC_neg_0006.d
QC_7	P2-H12	HILIC_neg.m	D:\MassHunter\Data\Demo\Demo_HILIC_neg\Demo_HILIC_neg_0007.d
QC_8	P2-H12	HILIC_neg.m	D:\MassHunter\Data\Demo\Demo_HILIC_neg\Demo_HILIC_neg_0008.d
QC_9	P2-H12	HILIC_neg.m	D:\MassHunter\Data\Demo\Demo_HILIC_neg\Demo_HILIC_neg_0009.d
QC_10	P2-H12	HILIC_neg.m	D:\MassHunter\Data\Demo\Demo_HILIC_neg\Demo_HILIC_neg_0010.d
A_1	P1-A1	HILIC_neg.m	D:\MassHunter\Data\Demo\Demo_HILIC_neg\Demo_HILIC_neg_0011.d
A_2	P1-A2	HILIC_neg.m	D:\MassHunter\Data\Demo\Demo_HILIC_neg\Demo_HILIC_neg_0012.d
A_3	P1-A3	HILIC_neg.m	D:\MassHunter\Data\Demo\Demo_HILIC_neg\Demo_HILIC_neg_0013.d
A_4	P1-A4	HILIC_neg.m	D:\MassHunter\Data\Demo\Demo_HILIC_neg\Demo_HILIC_neg_0014.d

The Sample Name column is identical to INTERNAL\_SAMPLE\_ID column in the processed sample information file. Sample position is generated automatically, and all QC positions are skipped. The Data File column includes the complete file paths. Wranglr will generate a subfolder for each mode.

A_165	P2-F9	HILIC_neg.m	D:\MassHunter\Data\Demo\Demo_HILIC_neg\Demo_HILIC_neg_0188.d
A_166	P2-F10	HILIC_neg.m	D:\MassHunter\Data\Demo\Demo_HILIC_neg\Demo_HILIC_neg_0189.d
A_67_auto_msms_10V	P1-F7	HILIC_neg_AutoMSMS_10V.m	D:\MassHunter\Data\Demo\Demo_HILIC_neg\Demo_HILIC_neg_0190.d
A_67_auto_msms_20V	P1-F7	HILIC_neg_AutoMSMS_20V.m	D:\MassHunter\Data\Demo\Demo_HILIC_neg\Demo_HILIC_neg_0191.d
A_67_auto_msms_40V	P1-F7	HILIC_neg_AutoMSMS_40V.m	D:\MassHunter\Data\Demo\Demo_HILIC_neg\Demo_HILIC_neg_0192.d
A_80_auto_msms_10V	P1-G8	HILIC_neg_AutoMSMS_10V.m	D:\MassHunter\Data\Demo\Demo_HILIC_neg\Demo_HILIC_neg_0193.d
A_80_auto_msms_20V	P1-G8	HILIC_neg_AutoMSMS_20V.m	D:\MassHunter\Data\Demo\Demo_HILIC_neg\Demo_HILIC_neg_0194.d
A_80_auto_msms_40V	P1-G8	HILIC_neg_AutoMSMS_40V.m	D:\MassHunter\Data\Demo\Demo_HILIC_neg\Demo_HILIC_neg_0195.d
QC_1_auto_msms_10V	P2-H12	HILIC_neg_AutoMSMS_10V.m	D:\MassHunter\Data\Demo\Demo_HILIC_neg\Demo_HILIC_neg_0196.d
QC_1_auto_msms_20V	P2-H12	HILIC_neg_AutoMSMS_20V.m	D:\MassHunter\Data\Demo\Demo_HILIC_neg\Demo_HILIC_neg_0197.d
QC_1_auto_msms_40V	P2-H12	HILIC_neg_AutoMSMS_40V.m	D:\MassHunter\Data\Demo\Demo_HILIC_neg\Demo_HILIC_neg_0198.d
STOP	Vial 1	STOP.m	D:\MassHunter\Data\Demo\Demo_HILIC_neg\Demo_HILIC_neg_0199.d

The AutoMSMS runs and the STOP run will be included in the end of the Worklist files. The position of STOP is always put to “Vial 1”

An individual worklist file will be generated for each mode. All the files generated by Wranglr can be downloaded as .csv files (with Windows line endings).

## MPP files tab

In the MPP tab you can combine MPP output files with sample information. The MPP output files should include the raw abundance values and the following columns: Compound, Mass, Retention Time, CompositeSpectrum, Frequency, CompoundAlgo and ionization mode.

The combined data matrix is ready for further data analysis. Other information of the compounds is stored in separate annotation file. Both files can be downloaded as .csv files.

Combine MPP output files

Project Title

Demo

Input processed sample info sheet

Browse...

sample\_info\_processed.csv

Upload complete

Upload the processed sample information file generated by Wranglr (.csv)

Sample information file OK

Found 664 original samples and 0 QC samples

The file is checked and the number of original and QC samples is calculated (this example file contains no QC samples)

Input HILIC NEG file

Browse...

Fold from MPP HILIC neg.txt

Upload complete

Upload the original MPP output files without any formatting

Input HILIC POS file

Browse...

Fold from MPP HILIC pos.txt

Upload complete

Input RP NEG file

Browse...

Fold from MPP RP neg.txt

Upload complete

Input RP POS file

Browse...

Fold from MPP RP pos.txt

Upload complete

Submit & Process

## Example of the combined data matrix:

Preview:

☒ Data matrix

☐ Annotation file

Show 25 entries

RUN_ORDER	INTERNAL_SAMPLE_ID	SAMPLE_ID	GROUP	QC	COMPOUND_835_5739a0_60039145_HILIC_neg	COMPOUND_880_6041a0_6022438_HILIC_neg
1	A_1	86	1	FALSE	36876	29488
2	A_2	59	2	FALSE	21472	26513
3	A_3	47	2	FALSE	37850	18807
4	A_4	12	2	FALSE	35907	26220
5	A_5	8	1	FALSE	39034	22130
6	A_6	60	1	FALSE	41822	27177
7	A_7	46	1	FALSE	31022	35585
8	A_8	25	3	FALSE	37226	31912
9	A_9	77	2	FALSE	64101	34478
10	A_10	100	3	FALSE	41231	24626
11	A_11	131	2	FALSE	42003	32209

The data matrix has one row per sample. The columns from the processed sample sheet are combined with abundance values of all compounds found in all modes. The compound names are modified to better suit data analysis pipeline.

Example of an annotation file:

Preview:

Data matrix

Annotation file

Show

25

entries

Search:

Compound	Mass	Retention.Time	CompositeSpectrum	Frequency	CompoundAlgo	Ionization.mode	chromatography
COMPOUND_835_5739a0_60039145_HILIC_neg	835.5739	0.6003915	(834.561, 618.15)(835.5636, 325.98)(836.5809, 2383.65)(837.5818, 1204.03)(838.5966, 6840.1)(839.6005, 3720.55)(840.5932, 2400.19)	161	FindByFormula	-	hilic
COMPOUND_880_6041a0_6022438_HILIC_neg	880.6041	0.6022438	(879.5964, 4105.43)(880.6052, 4976.04)(881.611, 2238.33)(882.6193, 1087.96)(883.6216, 516.88)	175	FindByFormula	-	hilic
COMPOUND_851_5713a0_61231047_HILIC_neg	851.5713	0.6123105	(850.5641, 40071.87)(851.5659, 20185.1)(852.5773, 24523.87)(853.5808, 11196.83)(854.5927, 18349.58)(855.5972, 8665.93)(856.6024, 4311.33)	181	FindByFormula	-	hilic
COMPOUND_827_5705a0_6225396_HILIC_neg	827.5705	0.6225396	(826.5634, 48164.44)(827.5662, 22517.98)(828.5764, 27366.94)(829.5799, 11596.97)(830.5916, 14789.17)(831.5959, 6514.86)(832.602, 3081.05)(833.6086, 1033.93)	181	FindByFormula	-	hilic
COMPOUND_1588_1497a0_63386375_HILIC_neg	1588.1497	0.6338637	(1587.1439, 1736.1)(1588.1519, 2181.12)(1589.1565, 1379.81)(1590.1635, 1079.31)(1591.1703, 599.48)	156	FindByFormula	-	hilic
COMPOUND_1583_1154a0_6332007_HILIC_neg	1583.1154	0.6332007	(1582.1046, 841.17)(1583.1107, 779.23)(1584.126, 2245.74)(1585.1256, 1989.31)(1586.1373, 2390.57)(1587.1411, 1810.44)(1588.1488, 1675.66)	165	FindByFormula	-	hilic
COMPOUND_803_5716a0_6422843_HILIC_neg	803.5716	0.6422843	(802.5643, 74265.74)(803.5665, 34034.41)(804.5776, 46775.65)(805.581, 19167.41)(806.5839, 5319.22)(807.5871, 1110.69)(808.5786, 257.74)	181	FindByFormula	-	hilic
COMPOUND_799_5405a0_6426328_HILIC_neg	799.5405	0.6426328	(798.5369, 502.62)(799.5325, 199.74)(800.5466, 7602.79)(801.5496, 3623.57)(1598.1224, 445.69)	179	FindByFormula	-	hilic
COMPOUND_743_5465a0_6734437_HILIC_neg	743.5465	0.6734437	(742.5393, 1332.7)(743.543, 587.2)(744.5544, 510.12)(745.5594, 150.68)	161	FindByFormula	-	hilic

The Chromatography column is added. Compound names correspond to those in the data matrix.